

NFκB-p65 (Acetyl-Lys218) Polyclonal Antibody

Catalog No: #HW122

Package Size: #HW122-1 50ul #HW122-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NFκB-p65 (Acetyl-Lys218) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-NFκB-p65 (K218) Polyclonal Antibody detects endogenous levels of NFκB-p65 protein only when acetylated at K218.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the human NFκB-p65 around the acetylation site of K218.
Target Name	NFκB-p65
Modification	Acetyl
Other Names	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Accession No.	Swiss-Prot: Q04206NCBI Gene ID: 5970
Uniprot	Q04206
GeneID	5970;
Target Species	human
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/20000

Not yet tested in other applications

Images



Western Blot analysis of MOUSE-BRAIN cells using Acetyl-NFκB-p65 (K218) Polyclonal Antibody



Western Blot analysis of JurkatK562 cells using Acetyl-NFκB-p65 (K218) Polyclonal Antibody

Note: This product is for in vitro research use only