

Recombinant Human Bcl-2-like protein 15(BCL2L15)

Catalog No: #AP76936



Package Size: #AP76936-1 20ug #AP76936-2 100ug #AP76936-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Bcl-2-like protein 15(BCL2L15)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-163aaSequence Info:Full Length of BC127719
Accession No.	Q5TBC7
Uniprot	Q5TBC7
GeneID	440603;
Calculated MW	44.7 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	MKSSQTFEEQTECIVNTLLMDFLSPTLQVASRNLCVDEVDVSGEPCSFDAIIAGRLRMLGDQFNNGELEASEAK NVIAETIKGQTGAILQNTVESLSKTWCAQDSSLAYERAFLAVSVKLELYMAHIAPEVVGQVAIPMTGMINGNQAI REFIQGQGGWENLES
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"Expression of pro-apoptotic Bcl-2 isoforms reduces during malignant transformation in the human gastrointestinal tract."

Dempsey C.E., Dive C., Fletcher D.J., Barnes F.A., Lobo A., Bingle C.D., Whyte M.K., Renshaw S.A.

FEBS Lett. 579:3646-3650(2005) Research Topic:Others

Note: This product is for in vitro research use only