

Histone H3 (Mono-Methyl-Lys79) Monoclonal Antibody

Catalog No: #HW157

Orders: order@signalwayantibody.com

Package Size: #HW157-1 50ul #HW157-2 100ul

Support: tech@signalwayantibody.com


Description

Product Name	Histone H3 (Mono-Methyl-Lys79) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Histone H3 (mono methyl K79) protein.
Immunogen Type	peptide
Immunogen Description	mono methyl peptide at K79 of Histone H3
Target Name	Histone H3
Modification	Methyl
Other Names	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/
Accession No.	Swiss-Prot: P68431/Q71DI3/P84243NCBI Gene ID: 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
Uniprot	P68431
GeneID	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-2000

Images



Western blot analysis of HeLa, diluted at

1) 1:1000

2) 1:3000

Note: This product is for in vitro research use only