

ECFP-Tag Monoclonal Antibody

Catalog No: #T616

Package Size: #T616-1 50ul #T616-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | ECFP-Tag Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Purification | The antibody was affinity-purified from antiserum by affinity-chromatography |
| Applications | WB |
| Specificity | The antibody detects endogenous ECFP-Tag protein. |
| Immunogen Type | Recombinant Protein |
| Immunogen Description | Recombinant Protein |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:2000-5000

Images



Western blot analysis of Recombinant ECFP Protein using ECFP-Tag Monoclonal Antibody.

Note: This product is for in vitro research use only