

Lamin B1 Monoclonal Antibody

Catalog No: #40413

Package Size: #40413-1 50ul #40413-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Lamin B1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IP
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous Lamin B1 protein.
Target Name	Lamin B1
Other Names	LMNB1; LMN2; LMNB; Lamin-B1
Accession No.	Swiss-Prot: P20700NCBI Gene ID: 4001
Uniprot	P20700
GeneID	4001;
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-5000

Immunoprecipitation:1:200


Images



1) Input: Mouse Brain Tissue Lysate

2) IP product: IP dilute 1:200

Western blot analysis of 1) HepG2, 2) 293T, 3) Mouse Brain Tissue, 4) Rat Brain Tissue, diluted at 1:5000.



Note: This product is for in vitro research use only