

## HSP70 Monoclonal Antibody

Catalog No: #40416

Package Size: #40416-1 50ul #40416-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	HSP70 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous HSP70 proteins.
Target Name	HSP70
Accession No.	Swiss-Prot: P34931/P08107NCBI Gene ID: 3305/3303/3304
Uniprot	P34931
GeneID	3305;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-2000

Immunofluorescence: 1:100-200

## Images



IF analysis of HeLa with antibody (Left) and DAPI (Right) diluted at 1:100.



Western blot analysis of 1) Hela, 2) Mouse Brain, diluted at 1:2000.

---

Note: This product is for in vitro research use only