

## CD45 Monoclonal Antibody

Catalog No: #40420

Package Size: #40420-1 50ul #40420-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CD45 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous CD45 proteins.
Target Name	CD45
Other Names	PTPRC; CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD45
Accession No.	Swiss-Prot: P08575NCBI Gene ID: 5788
Uniprot	P08575
GeneID	5788;
SDS-PAGE MW	147kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:2000

## Images



Western blot analysis of HeLa, diluted at 1:2000.

Note: This product is for in vitro research use only