

## XRCC4 Monoclonal Antibody

Catalog No: #40455

Package Size: #40455-1 50ul #40455-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	XRCC4 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IP
Species Reactivity	Hu
Specificity	The antibody detects endogenous XRCC4 proteins.
Target Name	XRCC4
Other Names	XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4
Accession No.	Swiss-Prot: Q13426NCBI Gene ID: 7518
Uniprot	Q13426
GeneID	7518;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:2000


Immunoprecipitation:1:200

## Images



1) Input: Hela Cell Lysate 2) IP product: IP dilute 1:200

Western blot analysis of 1) HeLa, 2) 293T, diluted at 1:3000.



Note: This product is for in vitro research use only