

## PR Monoclonal Antibody

Catalog No: #40481



Package Size: #40481-1 50ul #40481-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PR Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	IHC IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous PR protein.
Immunogen Type	peptide
Immunogen Description	Synthetic Peptide
Target Name	PR
Accession No.	Swiss-Prot: P06401NCBI Gene ID: 5241
Uniprot	P06401
GeneID	5241;
SDS-PAGE MW	38,43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:200 Immunofluorescence: 1:100-200

Note: This product is for in vitro research use only