

## Caspase-8 Monoclonal Antibody

Catalog No: #40483

Package Size: #40483-1 50ul #40483-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Caspase-8 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Caspase-8 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	Caspase-8
Accession No.	Swiss-Prot: Q14790NCBI Gene ID: 841
Uniprot	Q14790
GeneID	841;
SDS-PAGE MW	43,57kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-2000

Immunohistochemistry:1:200-500

## Images



Western blot analysis of 1) HeLa, 2) Mouse Brain Tissue, 3) Rat Brain Tissue using Caspase-8 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using Caspase-8 Monoclonal Antibody.

---

Note: This product is for in vitro research use only