

HAO1 Monoclonal Antibody

Catalog No: #40484

Package Size: #40484-1 50ul #40484-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HAO1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Ms Rt
Specificity	The antibody detects endogenous HAO1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	HAO1
Accession No.	Swiss-Prot: Q9UJM8NCBI Gene ID: 54363
Uniprot	Q9UJM8
GeneID	54363;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-2000

Images



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Monoclonal Antibody.

Note: This product is for in vitro research use only