

HAO1 Monoclonal Antibody

Catalog No: #40484

Package Size: #40484-1 50ul #40484-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | HAO1 Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Purification | The antibody was affinity-purified from antiserum by affinity-chromatography |
| Applications | WB |
| Species Reactivity | Ms Rt |
| Specificity | The antibody detects endogenous HAO1 protein. |
| Immunogen Type | Recombinant Protein |
| Immunogen Description | Recombinant Protein |
| Target Name | HAO1 |
| Accession No. | Swiss-Prot: Q9UJM8NCBI Gene ID: 54363 |
| Uniprot | Q9UJM8 |
| GeneID | 54363; |
| SDS-PAGE MW | 41kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:1000-2000

Images



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Monoclonal Antibody.

Note: This product is for in vitro research use only