

Survivin Monoclonal Antibody

Catalog No: #40488

Package Size: #40488-1 50ul #40488-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Survivin Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous Survivin protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	Survivin
Accession No.	Swiss-Prot: O15392NCBI Gene ID: 332
Uniprot	O15392
GeneID	332;
SDS-PAGE MW	16kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000-2000

Immunohistochemistry: 1:200-500

Images



Western blot analysis of 1) HeLa, 2) 293, 3) PC12 using Survivin Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma using Survivin Monoclonal Antibody.

Note: This product is for in vitro research use only