

CREB-1 Monoclonal Antibody

Catalog No: #40489

Package Size: #40489-1 50ul #40489-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CREB-1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous CREB-1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	CREB-1
Accession No.	Swiss-Prot: P16220NCBI Gene ID: 1385
Uniprot	P16220
GeneID	1385;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:500-1000

Immunohistochemistry: 1:100-200

Images



Western blot analysis of 1) A431, 2) 3T3 using CREB-1 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using CREB-1 Monoclonal Antibody.

Note: This product is for in vitro research use only