

MICU1 Monoclonal Antibody

Catalog No: #40490

Package Size: #40490-1 50ul #40490-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MICU1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous MICU1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	MICU1
Accession No.	Swiss-Prot: Q9BPX6NCBI Gene ID: 10367
Uniprot	Q9BPX6
GeneID	10367;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:1000-2000

Immunohistochemistry: 1:100-200

Images



Western blot analysis of 1) MCF7, 2) Mouse Brain Tissue, 3) Rat Brain Tissue using MICU1 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using MICU1 Monoclonal Antibody.

Note: This product is for in vitro research use only