

Active Caspase-3 Monoclonal Antibody

Catalog No: #40491

Package Size: #40491-1 50ul #40491-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Active Caspase-3 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Active Caspase-3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	Active Caspase-3
Accession No.	Swiss-Prot: P42574NCBI Gene ID: 836
Uniprot	P42574
GeneID	836;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:500-1000

Immunohistochemistry: 1:100-200

Images



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Brain Tissue using Active Caspase-3 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Active Caspase-3 Monoclonal Antibody.

Note: This product is for in vitro research use only