

GAPDH Monoclonal Antibody

Catalog No: #40493

Package Size: #40493-1 50ul #40493-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GAPDH Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu Rt Ms Mk Dg C Hm Rb Pg Sh
Specificity	The antibody detects endogenous GAPDH protein.
Target Name	GAPDH
Other Names	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Accession No.	Swiss-Prot: P04406NCBI Gene ID: 2597
Uniprot	P04406
GeneID	2597;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-5000

Images



Western blot analysis of HeLa, diluted at 1:3000.

Note: This product is for in vitro research use only