

α -SMA Monoclonal Antibody

Catalog No: #40495

Package Size: #40495-1 50ul #40495-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	α -SMA Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous α -SMA protein.
Immunogen Type	peptide
Immunogen Description	Synthetic Peptide
Target Name	α -SMA
Accession No.	Swiss-Prot: P62736NCBI Gene ID: 59
Uniprot	P62736
GeneID	59;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:10000-100000

Immunohistochemistry:1:1000-2000

Images



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Brain using α -SMA Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Mouse Cecal Tissue using α -SMA Monoclonal Antibody.

Note: This product is for in vitro research use only