

Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody

Catalog No: #40504

Package Size: #40504-1 50ul #40504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-9 p35 protein resulting from cleavage adjacent to D315.
Immunogen Description	Synthesized peptide derived from the Internal region of human Caspase-9 p35.
Target Name	Cleaved-Caspase-9 p35
Other Names	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6
Accession No.	Swiss-Prot: P55211NCBI Gene ID: 842
Uniprot	P55211
GeneID	842;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody



Western Blot analysis of VEC A549 cells using Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody

Note: This product is for in vitro research use only