

Adiponectin Polyclonal Antibody

Catalog No: #40505

Package Size: #40505-1 50ul #40505-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adiponectin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Adiponectin Polyclonal Antibody detects endogenous levels of Adiponectin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Adiponectin.
Target Name	Adiponectin
Other Names	ADIPOQ; ACDC; ACRP30; APM1; GBP28; Adiponectin; 30 kDa adipocyte complement-related protein; Adipocyte complement-related 30 kDa protein; ACRP30; Adipocyte; C1q and collagen domain-containing protein; Adipose most abundant gene transcript 1
Accession No.	Swiss-Prot: Q15848NCBI Gene ID: 9370
Uniprot	Q15848
GeneID	9370;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of A549 cells using Adiponectin Polyclonal Antibody



Note: This product is for in vitro research use only