

ABCA8 Polyclonal Antibody

Catalog No: #40532

Package Size: #40532-1 50ul #40532-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABCA8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	ABCA8 Polyclonal Antibody detects endogenous levels of ABCA8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ABCA8.
Target Name	ABCA8
Other Names	ABCA8; KIAA0822; ATP-binding cassette sub-family A member 8
Accession No.	Swiss-Prot: O94911NCBI Gene ID: 10351
Uniprot	O94911
GeneID	10351;
SDS-PAGE MW	179kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HuvEc cells using ABCA8 Polyclonal Antibody



Western Blot analysis of HuvEc PC12 cells using ABCA8 Polyclonal Antibody

Note: This product is for in vitro research use only