

ABCC12 Polyclonal Antibody

Catalog No: #40533

Package Size: #40533-1 50ul #40533-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCC12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ABCC12 Polyclonal Antibody detects endogenous levels of ABCC12 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ABCC12.
Target Name	ABCC12
Other Names	ABCC12; MRP9; Multidrug resistance-associated protein 9; ATP-binding cassette sub-family C member 12
Accession No.	Swiss-Prot: Q96J65NCBI Gene ID: 94160
Uniprot	Q96J65
GeneID	94160;
SDS-PAGE MW	170kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using ABCC12 Polyclonal Antibody

Note: This product is for in vitro research use only