

## ABCD1 Polyclonal Antibody

Catalog No: #40535

Package Size: #40535-1 50ul #40535-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ABCD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ABCD1 Polyclonal Antibody detects endogenous levels of ABCD1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ABCD1.
Target Name	ABCD1
Other Names	ABCD1; ALD; ATP-binding cassette sub-family D member 1; Adrenoleukodystrophy protein; ALDP
Accession No.	Swiss-Prot: P33897NCBI Gene ID: 215
Uniprot	P33897
GeneID	215;
SDS-PAGE MW	75kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of HeLa cells using ABCD1 Polyclonal Antibody

---

Note: This product is for in vitro research use only