

ABCD4 Polyclonal Antibody

Catalog No: #40536

Package Size: #40536-1 50ul #40536-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCD4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	ABCD4 Polyclonal Antibody detects endogenous levels of ABCD4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ABCD4.
Target Name	ABCD4
Other Names	ABCD4; PXMP1L; ATP-binding cassette sub-family D member 4; PMP70-related protein; P70R; Peroxisomal membrane protein 1-like; PXMP1-L; Peroxisomal membrane protein 69; PMP69
Accession No.	Swiss-Prot: O14678NCBI Gene ID: 5826
Uniprot	O14678
GeneID	5826;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using ABCD4 Polyclonal Antibody



Western Blot analysis of K562 cells using ABCD4 Polyclonal Antibody

Note: This product is for in vitro research use only