

Abi-1 Polyclonal Antibody

Catalog No: #40538

Package Size: #40538-1 50ul #40538-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Abi-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Abi-1 Polyclonal Antibody detects endogenous levels of Abi-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Abi-1.
Target Name	Abi-1
Other Names	ABI1; SSH3BP1; Abl interactor 1; Abelson interactor 1; Abi-1; Abl-binding protein 4; AblBP4; Eps8 SH3 domain-binding protein; Eps8-binding protein; Nap1-binding protein; Nap1BP; Spectrin SH3 domain-binding protein 1; e3B1
Accession No.	Swiss-Prot: Q8IZP0NCBI Gene ID: 10006
Uniprot	Q8IZP0
GeneID	10006;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Abi-1 Polyclonal Antibody



Western Blot analysis of HuvEc cells using Abi-1 Polyclonal Antibody

Note: This product is for in vitro research use only