

AChE Polyclonal Antibody

Catalog No: #40541

Package Size: #40541-1 50ul #40541-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AChE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	AChE Polyclonal Antibody detects endogenous levels of AChE protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human AChE.
Target Name	AChE
Other Names	ACHE; Acetylcholinesterase; AChE
Accession No.	Swiss-Prot: P22303NCBI Gene ID: 43
Uniprot	P22303
GeneID	43;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Immunohistochemistry: 1/100 - 1/300

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using AChE Polyclonal Antibody

Note: This product is for in vitro research use only