

ACOT2 Polyclonal Antibody

Catalog No: #40543

Package Size: #40543-1 50ul #40543-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACOT2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu
Specificity	ACOT2 Polyclonal Antibody detects endogenous levels of ACOT2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ACOT2.
Target Name	ACOT2
Other Names	ACOT2; PTE2; PTE2A; Acyl-coenzyme A thioesterase 2; mitochondrial; Acyl-CoA thioesterase 2; Acyl-coenzyme A thioester hydrolase 2a; CTE-1a; Long-chain acyl-CoA thioesterase 2; ZAP128
Accession No.	Swiss-Prot: P49753NCBI Gene ID: 10965
Uniprot	P49753
GeneID	10965;
SDS-PAGE MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of VEC cells using ACOT2 Polyclonal Antibody

Note: This product is for in vitro research use only