

ACSL6 Polyclonal Antibody

Catalog No: #40547

Package Size: #40547-1 50ul #40547-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACSL6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	ACSL6 Polyclonal Antibody detects endogenous levels of ACSL6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ACSL6.
Target Name	ACSL6
Other Names	ACSL6; ACS2; FACL6; KIAA0837; LACS5; Long-chain-fatty-acid--CoA ligase 6; Long-chain acyl-CoA synthetase 6; LACS 6
Accession No.	Swiss-Prot: Q9UKU0NCBI Gene ID: 23305
Uniprot	Q9UKU0
GeneID	23305;
SDS-PAGE MW	78kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using ACSL6 Polyclonal Antibody

Note: This product is for in vitro research use only