

Actin β Polyclonal Antibody

Catalog No: #40552

Package Size: #40552-1 50ul #40552-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Actin β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Actin β Polyclonal Antibody detects endogenous levels of Actin β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Actin β .
Target Name	Actin β
Other Names	ACTB; Actin; cytoplasmic 1; Beta-actin
Accession No.	Swiss-Prot: P60709NCBI Gene ID: 60
Uniprot	P60709
GeneID	60;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/1000 - 1/4000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of HepG2 cells using Actin β Polyclonal Antibody

Western Blot analysis of HELA cells using Actin β Polyclonal Antibody



Note: This product is for in vitro research use only