

Adducin α/β Polyclonal Antibody

Catalog No: #40554

Package Size: #40554-1 50ul #40554-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adducin α/β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IP IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Adducin α/β Polyclonal Antibody detects endogenous levels of Adducin α/β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Adducin α/β around the non-phosphorylation site of S726/713.
Target Name	Adducin α/β
Other Names	ADD1; ADDA; Alpha-adducin; Erythrocyte adducin subunit alpha; ADD2; ADDB; Beta-adducin; Erythrocyte adducin subunit beta
Accession No.	Swiss-Prot: P35611/P35612NCBI Gene ID: 118/119
Uniprot	P35611
GeneID	118;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunoprec


Immunoprecipitation: 2-5 ug/mg lysate.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HELA cells using Adducin α/β Polyclonal Antibody



Western Blot analysis of HELA PC12 NIH-3T3 cells using Adducin α/β Polyclonal Antibody

Note: This product is for in vitro research use only