

Adenosine A3-R Polyclonal Antibody

Catalog No: #40555

Package Size: #40555-1 50ul #40555-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine A3-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Adenosine A3-R Polyclonal Antibody detects endogenous levels of Adenosine A3-R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Adenosine A3-R.
Target Name	Adenosine A3-R
Other Names	ADORA3; Adenosine receptor A3
Accession No.	Swiss-Prot: P33765NCBI Gene ID: 140
Uniprot	P0DMS8
GeneID	140;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using Adenosine A3-R Polyclonal Antibody



Western Blot analysis of MCF7 HELA PC12 cells using Adenosine A3-R Polyclonal Antibody

Note: This product is for in vitro research use only