

AIM2 Polyclonal Antibody

Catalog No: #40561

Package Size: #40561-1 50ul #40561-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AIM2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	AIM2 Polyclonal Antibody detects endogenous levels of AIM2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human AIM2.
Target Name	AIM2
Other Names	AIM2; Interferon-inducible protein AIM2; Absent in melanoma 2
Accession No.	Swiss-Prot: O14862NCBI Gene ID: 9447
Uniprot	O14862
GeneID	9447;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/40000.

Not yet tested in other applications.


Images




Western Blot analysis of HuvEc cells using AIM2 Polyclonal Antibody



Western Blot analysis of HuvEc cells using AIM2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human tonsil.



Immunofluorescence analysis of Hela cell.

Note: This product is for in vitro research use only