

AIRE-1 Polyclonal Antibody

Catalog No: #40562

Package Size: #40562-1 50ul #40562-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AIRE-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms
Specificity	AIRE-1 Polyclonal Antibody detects endogenous levels of AIRE-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human AIRE-1.
Target Name	AIRE-1
Other Names	AIRE; APECED; Autoimmune regulator; Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein; APECED protein
Accession No.	Swiss-Prot: O43918NCBI Gene ID: 326
Uniprot	O43918
GeneID	326;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of 293 cells using AIRE-1 Polyclonal Antibody



Note: This product is for in vitro research use only