

AKAP 95 Polyclonal Antibody

Catalog No: #40564

Package Size: #40564-1 50ul #40564-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AKAP 95 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	AKAP 95 Polyclonal Antibody detects endogenous levels of AKAP 95 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human AKAP 95.
Target Name	AKAP 95
Other Names	AKAP8; AKAP95; A-kinase anchor protein 8; AKAP-8; A-kinase anchor protein 95 kDa; AKAP 95
Accession No.	Swiss-Prot: O43823NCBI Gene ID: 10270
Uniprot	O43823
GeneID	10270;
SDS-PAGE MW	76kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using AKAP 95 Polyclonal Antibody



Western Blot analysis of PC12 MCF7 Jurkat cells using AKAP 95 Polyclonal Antibody

Note: This product is for in vitro research use only