

ALDH1B1 Polyclonal Antibody

Catalog No: #40579

Package Size: #40579-1 50ul #40579-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ALDH1B1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Mk
Specificity	ALDH1B1 Polyclonal Antibody detects endogenous levels of ALDH1B1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ALDH1B1.
Target Name	ALDH1B1
Other Names	ALDH1B1; ALDH5; ALDHX; Aldehyde dehydrogenase X; mitochondrial; Aldehyde dehydrogenase 5; Aldehyde dehydrogenase family 1 member B1
Accession No.	Swiss-Prot: P30837NCBI Gene ID: 219
Uniprot	P30837
GeneID	219;
SDS-PAGE MW	57kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using ALDH1B1 Polyclonal Antibody



Western Blot analysis of HELA cells using ALDH1B1 Polyclonal Antibody

Note: This product is for in vitro research use only