

ALDH3B1 Polyclonal Antibody

Catalog No: #40580

Package Size: #40580-1 50ul #40580-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ALDH3B1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	ALDH3B1 Polyclonal Antibody detects endogenous levels of ALDH3B1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ALDH3B1.
Target Name	ALDH3B1
Other Names	ALDH3B1; ALDH7; Aldehyde dehydrogenase family 3 member B1; Aldehyde dehydrogenase 7
Accession No.	Swiss-Prot: P43353NCBI Gene ID: 221
Uniprot	P43353
GeneID	221;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using ALDH3B1 Polyclonal Antibody

Note: This product is for in vitro research use only