

## Aldolase B Polyclonal Antibody

Catalog No: #40582

Package Size: #40582-1 50ul #40582-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Aldolase B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Aldolase B Polyclonal Antibody detects endogenous levels of Aldolase B protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Aldolase B.
Target Name	Aldolase B
Other Names	ALDOB; ALDB; Fructose-bisphosphate aldolase B; Liver-type aldolase
Accession No.	Swiss-Prot: P05062NCBI Gene ID: 229
Uniprot	P05062
GeneID	229;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HuvEc cells using Aldolase B Polyclonal Antibody



Western Blot analysis of HuvEc MCF7 HELA PC12 cells using Aldolase B Polyclonal Antibody

Note: This product is for in vitro research use only