

AOX1 Polyclonal Antibody

Catalog No: #40597

Package Size: #40597-1 50ul #40597-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AOX1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Rt
Specificity	AOX1 Polyclonal Antibody detects endogenous levels of AOX1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human AOX1.
Target Name	AOX1
Other Names	AOX1; AO; Aldehyde oxidase
Accession No.	Swiss-Prot: Q06278NCBI Gene ID: 316
Uniprot	Q06278
GeneID	316;
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using AOX1 Polyclonal Antibody

Note: This product is for in vitro research use only