

## AP-1 Polyclonal Antibody

Catalog No: #40598

Package Size: #40598-1 50ul #40598-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AP-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IP ELISA
Species Reactivity	Hu Ms Rt
Specificity	AP-1 Polyclonal Antibody detects endogenous levels of AP-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human AP-1 around the non-phosphorylation site of Y170.
Target Name	AP-1
Other Names	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
Accession No.	Swiss-Prot: P05412NCBI Gene ID: 3725
Uniprot	P05412
GeneID	3725;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300. Immunoprec

Immunoprecipitation: 2-5 ug/mg lysate.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293 cells using AP-1 Polyclonal Antibody

---

Note: This product is for in vitro research use only