

APLP-2 Polyclonal Antibody

Catalog No: #40603

Package Size: #40603-1 50ul #40603-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	APLP-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	APLP-2 Polyclonal Antibody detects endogenous levels of APLP-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human APLP-2.
Target Name	APLP-2
Other Names	APLP2; APPL2; Amyloid-like protein 2; APLP-2; APPH; Amyloid protein homolog; CDEI box-binding protein; CDEBP
Accession No.	Swiss-Prot: Q06481NCBI Gene ID: 334
Uniprot	Q06481
GeneID	334;
SDS-PAGE MW	87kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using APLP-2 Polyclonal Antibody



Western Blot analysis of 293 NIH-3T3 cells using APLP-2 Polyclonal Antibody

Note: This product is for in vitro research use only