

ApoE Polyclonal Antibody

Catalog No: #40604

Package Size: #40604-1 50ul #40604-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ApoE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	ApoE Polyclonal Antibody detects endogenous levels of ApoE protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human ApoE.
Target Name	ApoE
Other Names	APOE; Apolipoprotein E; Apo-E
Accession No.	Swiss-Prot: P02649NCBI Gene ID: 348
Uniprot	P02649
GeneID	348;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using ApoE Polyclonal Antibody

Note: This product is for in vitro research use only