

ApoF Polyclonal Antibody

Catalog No: #40605

Package Size: #40605-1 50ul #40605-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ApoF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ApoF Polyclonal Antibody detects endogenous levels of ApoF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human ApoF.
Target Name	ApoF
Other Names	APOF; Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP
Accession No.	Swiss-Prot: Q13790NCBI Gene ID: 319
Uniprot	Q13790
GeneID	319;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using ApoF Polyclonal Antibody



Western Blot analysis of KB PC-3 MCF7 cells using ApoF Polyclonal Antibody

Note: This product is for in vitro research use only