

AR- $\beta$ 3 Polyclonal Antibody

Catalog No: #40617

Package Size: #40617-1 50ul #40617-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AR- $\beta$ 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	AR- $\beta$ 3 Polyclonal Antibody detects endogenous levels of AR- $\beta$ 3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human AR- $\beta$ 3.
Target Name	AR- $\beta$ 3
Other Names	ADRB3; ADRB3R; B3AR; Beta-3 adrenergic receptor; Beta-3 adrenoreceptor; Beta-3 adrenoceptor
Accession No.	Swiss-Prot: P13945NCBI Gene ID: 155
Uniprot	P13945
GeneID	155;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

## Images

Western Blot analysis of K562 cells using AR- $\beta$ 3 Polyclonal Antibody

Western Blot analysis of 22RV1 cells using AR- $\beta$ 3 Polyclonal Antibody



Note: This product is for in vitro research use only