

Bone morphogenetic protein 4

Catalog No: #AP78322



Package Size: #AP78322-1 50ug #AP78322-2 100ug #AP78322-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone morphogenetic protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	20-292AA
Other Names	Bmp-4,Dvr-4,Bone morphogenetic protein 2B,BMP-2B
Accession No.	P21275Gene name:Bmp4
Uniprot	P21275
GeneID	12159;
Calculated MW	30.03
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Induces cartilage and bone formation. Acts in concert with PTHLH/PTHRP to stimulate ductal outgrowth during embryonic mammary development and to inhibit hair follicle induction.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only