

Interleukin-12 subunit alpha

Catalog No: #AP78347



Package Size: #AP78347-1 50ug #AP78347-2 100ug #AP78347-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-12 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	26-221AA
Other Names	Cytotoxic lymphocyte maturation factor 35 kDa subunit,IL-12 subunit p35,CLMF p35
Accession No.	P54349Gene name:IL12A
Uniprot	P54349
GeneID	281856;
Calculated MW	21.56
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only