

Interleukin-18

Catalog No: #AP78351

Package Size: #AP78351-1 50ug #AP78351-2 100ug #AP78351-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	37-193AA
Other Names	IGIF,IL1F4,Ibctadekin,Interferon gamma-inducing factor,Interleukin-1 gamma,IFN-gamma-inducing factor,IL-1 gamma
Accession No.	Q14116Gene name:IL18
Uniprot	Q14116
GeneID	3606;
Calculated MW	17.27
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Augments natural killer cell activity in spleen cells and stimulates interferon gamma production in T-helper type I cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only