

Ataxin-1 Polyclonal Antibody

Catalog No: #40621

Package Size: #40621-1 50ul #40621-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ataxin-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	Ataxin-1 Polyclonal Antibody detects endogenous levels of Ataxin-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Ataxin-1 around the non-phosphorylation site of S776.
Target Name	Ataxin-1
Other Names	ATXN1; ATX1; SCA1; Ataxin-1; Spinocerebellar ataxia type 1 protein
Accession No.	Swiss-Prot: P54253NCBI Gene ID: 6310
Uniprot	P54253
GeneID	6310;
SDS-PAGE MW	87kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Images



Western Blot analysis of MOUSE-BRAIN cells using Ataxin-1 Polyclonal Antibody

Note: This product is for in vitro research use only