

Interleukin-5

Catalog No: #AP78361



Package Size: #AP78361-1 50ug #AP78361-2 100ug #AP78361-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	20-132AA
Other Names	IL-5,B-cell growth factor II,Cytotoxic T-lymphocyte inducer,Eosinophil differentiation factor,T-cell replacing factor,BCGF-II,TRF
Accession No.	Q08125Gene name:IL5
Uniprot	Q08125
Calculated MW	12.43
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Factor that induces terminal differentiation of late-developing B-cells to immunoglobulin secreting cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only