

## Leptin receptor

Catalog No: #AP78369

Package Size: #AP78369-1 50ug #AP78369-2 100ug #AP78369-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leptin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	863-1165AA
Other Names	DB,OBR,HuB219,OB receptor,OB-R,CD295
Accession No.	P48357Gene name:LEPR
Uniprot	P48357
GeneID	3953;
Calculated MW	33.33
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Isoform E: Antagonizes Isoform A and isoform B-mediated LEP binding and endocytosis.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only