

Ataxin-2L Polyclonal Antibody

Catalog No: #40623

Package Size: #40623-1 50ul #40623-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Ataxin-2L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	Ataxin-2L Polyclonal Antibody detects endogenous levels of Ataxin-2L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ataxin-2L.
Target Name	Ataxin-2L
Other Names	ATXN2L; A2D; A2LG; A2LP; A2RP; Ataxin-2-like protein; Ataxin-2 domain protein; Ataxin-2-related protein
Accession No.	Swiss-Prot: Q8WWM7NCBI Gene ID: 11273
Uniprot	Q8WWM7
GeneID	11273;
SDS-PAGE MW	113kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 293T cells using Ataxin-2L Polyclonal Antibody

Note: This product is for in vitro research use only